waterpik

Congratulations!

You are now on your way to better oral health!

WE'RE HERE TO HELP

We hope you are loving your new Toothbrush. If for any reason you are not satisfied or have any product questions, please contact our fast and friendly Customer Service Team based in Fort Collins, Colorado.

Call: 1-800-525-2774
Monday-Friday
7 am - 5 pm MST
Chat or E-mail
service@waterpik.com

INSTRUCTION MANUAL INSIDE
FLIP OVER FOR QUICK START GUIDE

IMPORTANT SAFEGUARDS

READ ALL INSTRUCTIONS BEFORE USING.

- This product is for household use.

When using electrical products, especially when children are present, basic safety precautions should always be followed, including the following:

DANGER:

To reduce the risk of electrocution:

- Do not use while bathing.
Do not place or store product where it can fall or be pulled into a tub or sink.
Do not place or drop the product into water or other liquid.
Do not reach for a product that has fallen into water. Unplug immediately.
Do not handle the user supplied USB power supply (power adapter) with wet hands.

- Unplug the user supplied USB power supply (power adapter) when not in use.
Check the charging cord and user supplied USB power supply (power adapter) for damage before the first use and during the life of the product.
Do not use this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into any liquid. Contact Water Pik, Inc. at 1-800-525-2774 (US) 1-888-226-3042 (Canada) or visit www.waterpik.com (US) www.waterpik.ca (Canada).

WARNING:

To reduce the risk of burns, electrocution, fire, or injury to persons:

- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
Children should be supervised to ensure that they do not play with the appliance.
Use this product only for its intended use as described in these instructions or as recommended by your dental professional.
Only use replacement brush heads and accessories that are recommended by Water Pik, Inc. as identified in this instruction manual.
If the charging cord is damaged, it must be replaced by a special cord or assembly available from the manufacturer or its service agent.

- If your physician advised you to receive antibiotic premedication before dental procedures, you should consult your dentist before using this instrument or any other oral hygiene aid.
Remove any oral jewelry prior to use of this product.
Do not use if you have an open wound on your tongue or in your mouth.
There are no consumer serviceable electrical parts in this product and it does not require maintenance.
Ensure contacts on the cord, the unit and the power adapter are clean and dry while the unit is being charged.
Do not expose the battery or appliance to excessive temperatures.
Using toothpastes containing peroxide may damage the plastic housing of this product.

Limited Two Year Warranty

Water Pik, Inc. warrants to the original consumer of this new product that it is free from defects in materials and workmanship for 2 years from the date of purchase. Consumer will be required to submit the original purchase receipt as proof of purchase date and if requested, the entire product, to support a warranty claim. Water Pik, Inc. will replace at its discretion the product, provided the product has not been abused, misused, altered or damaged after purchase, was used according to instructions, and was used only with accessories or consumable parts approved by Water Pik, Inc. This limited warranty excludes accessories or consumable parts such as tips, etc. Use with improper voltage will damage product and VOID warranty.

This warranty gives you specific legal rights. You may also have other rights that vary, depending on the law in your place of residence.

WATER PIK, INC. DISCLAIMS ALL OTHER EXPRESS OR IMPLIED WARRANTIES TO THE EXTENT PERMITTED BY LAW INCLUDING, BUT NOT LIMITED TO, ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. ALL IMPLIED WARRANTIES ARE LIMITED TO THE DURATION OF THIS LIMITED WARRANTY. SOME STATES OR PROVINCES MAY NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, SO THE ABOVE MAY NOT APPLY TO YOU.

SAVE THESE INSTRUCTIONS

Have Questions or Need More Information?

- FULL INSTRUCTION MANUAL
CLEANING
REPLACEMENT BRUSH HEADS
TROUBLESHOOTING
COMMON QUESTIONS

Let Us Help You Get Started!

We have many helpful tools to provide you with the best possible experience. Visit: waterpik.com/sensonic (USA) waterpik.ca/sensonic (Canada)

SERVICE MAINTENANCE

Waterpik™ electric toothbrushes have no consumer serviceable electrical items and do not require routine service maintenance. For all parts and accessories go to www.waterpik.com (US) or www.waterpik.ca (Canada). You may also call our toll-free customer service number at 1-800-525-2774 (US) or 1-888-226-3042 (Canada).

Refer to the serial and model numbers in all correspondence. These numbers are located on the base or handle of your toothbrush.

Still have questions? Go to www.waterpik.com (US) or www.waterpik.ca (Canada)

CLEANING

- 1. After each use, turn off unit, remove the brush head from the handle and rinse brush head under running water.
2. Wipe the handle dry after each use.
3. Periodically, clean the charger base, first unplugging the unit and then wiping with a damp cloth.

BATTERY CARE

If the Water Flosser is to be stored for an extended period of time (six months or more), be sure to fully charge it prior to storage.

Do not expose the battery or appliance to excessive temperatures.

Limitation of Liability

This limited warranty is your sole remedy against Water Pik, Inc. and it shall not be liable for any damages, whether direct, indirect, incidental, special, consequential, exemplary or otherwise, including but not limited to property damage, arising out of any theory of recovery, including statutory, contract or tort.

Arbitration

BINDING ARBITRATION AND WAIVER OF RIGHT TO ASSERT MATTERS IN COURT AND RIGHT TO JURY TRIAL.

By purchasing, installing, or using this product you agree that any controversy or claim arising out of or in relation to your purchase, use, modification, installation, or removal of this product shall be settled by binding arbitration. The arbitrator's decision will be final and binding. Judgment on the award rendered by the arbitrator may be entered in any court having jurisdiction over the matter. Further, you agree that you will not have the right to litigate any claim arising out of or in relation to your purchase, use, modification, installation, or removal of this product in any court or to have a jury trial on that claim.

Arbitration of any matters shall be administered by the American Arbitration Association and conducted by one arbitrator in accordance with the consumer arbitration rules of the American Arbitration Association in effect at the time of filing of the demand. The venue for all such arbitrations shall be Colorado

and any such arbitration shall be subject to the Federal Arbitration Act and Colorado law.

WAIVER OF RIGHT TO PARTICIPATE IN CLASS ACTIONS OR PURSUE CLAIMS IN A REPRESENTATIVE CAPACITY.

Further, by purchasing, installing, or using this product you agree that arbitration must be on an individual basis. This means neither you, nor we, may join or consolidate claims in arbitration by or against other purchasers or users of Water Pik, Inc. products or the general public.

Additionally, you may not litigate in court or arbitrate any claims as a representative or member of a class or in a purported representative capacity on behalf of the general public, other purchasers or users of this product, other persons or entities similarly situated, or in a private attorney general capacity.

FCC Statement (USA)

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1. This device may not cause harmful interference, and
2. This device must accept any interference received, including interference that may cause undesired operation.

Warning: Changes or modifications not expressly approved by Water Pik, Inc. could void the user's authority to operate the equipment.

- 3 Remove and discard the brush head. Cover the handle entirely with the towel or cloth.
4 Clamp the top of the handle against your work surface with one hand and strike the handle housing around the bottom perimeter on all sides to release the end cap latch. Continue striking and rotating the handle until the end cap releases.
5 Remove the internal components by turning the handle upside down, then un-latch the plastic tabs and slide the chassis from the outer housing.
6 The internal components should slide out when the handle is positioned upside-down. If this does not occur, repeat steps 5 and 6 until the components release.
7 Remove the battery from the chassis. Cut the battery wires one at a time with the wire cutters.
8 Recycle and dispose of the battery properly. The remaining product pieces should be discarded appropriately.

WARNING:

DO NOT cut both battery wires at the same time. DO NOT disassemble, puncture, hit, squeeze or throw the battery into fire.

- 1 Run down the battery by operating the BRUSH button until the motor does not run.
2 Obtain a towel or cloth, hammer, flat-head screwdriver, and wire cutters.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

IC Statement (Canada)



This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

- 1. This device may not cause interference.
2. This device must accept any interference, including interference that may cause undesired operation of the device.

Quick Start Guide

GETTING STARTED


- 1 Charge for 4 hours before first use. Light ring will flash on handle while charging. When fully charged all lights will shut off. Light ring will glow red when battery is running low.

- 2 Slide brush head onto handle and push it down firmly.


- 3 Place brush head into your mouth before turning unit on.

- 5 For best results, position brush along gumline at a 45 degree angle. Use light pressure with a slow back and forth motion covering one or two teeth at a time.

- 6 Bristles change from green to yellow when time to replace. (For optimal performance replace brush heads every 3 months.)

Replacement Brush Heads:

Sensonic™ Brush Head – Replace every 3 months

Replacement Brush Heads may be purchased online at www.waterpik.com (US) www.waterpik.ca (Canada), or by phone at 1-800-525-2774 (US) 1-888-226-3042 (Canada).

Rated Input: 5.0 VDC, 1.0 A Max
Rechargeable Battery: Lithium-ion

Table with 7 columns: Symbol, Definition, Symbol, Definition, Symbol, Definition, Symbol, Definition. It lists various safety and regulatory symbols like 'Refer to instruction manual', 'Detachable Power Supply', 'Manufacturer', 'Direct Current (DC) Voltage', 'Ingress Protection Code', 'IEC Class II Equipment', 'Do not use this device in a shower, bathtub and do not immerse in any liquid', 'Authorized Representative in the European Community', 'CE marking of conformity', 'UK Conformity Assessed', and 'Regulatory Compliance Mark'.

Call: 1-800-525-2774 (USA) or 1-888-226-3042 (Canada)
Monday-Friday: 7 am - 5 pm MST
Chat or E-mail at service@waterpik.com

CONNECT WITH US @WATERPIK

Call: 1-800-525-2774 (USA) or 1-888-226-3042 (Canada)
www.waterpik.com or www.waterpik.ca

waterpik

Water Pik, Inc., a subsidiary of Church & Dwight Co., Inc.
1730 East Prospect Road
Fort Collins, CO 80553-0001 USA
1-800-525-2774 (USA/E.-U.)
1-888-226-3042 (Canada)
www.waterpik.com
www.waterpik.ca
Designed and Tested in the USA.
Form No. 20033027-F AA
©2024 Water Pik, Inc.

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Conçu et testé aux États-Unis. Numéro du modèle 20033027-FAA ©2024 Water Pik, Inc.