


Quick Start Guide



For use with all Waterpik® Water Flosser Models except: Sidekick (WF-04), Classic Professional (WP-72), Sonic-Fusion (SF-01/SF-02), Sonic-Fusion 2.0 (SF-03/SF-04), and Whitening Water Flossers (WF-05/WF-06)

NOTE: When used with handheld cordless Waterpik® Water Flossers, the tablet may not completely dissolve if only using 1 reservoir. If this happens, do one of the following:

- Extend water flossing time
- Discard excess tablet
- Leave leftover tablet in tip



We're Here to Help

We hope you are loving your new Boost Tip + Whitening Tablets.

If for any reason you are not satisfied or have any product questions, please contact our fast and friendly Customer Service Team based in Fort Collins, Colorado.



Call: 1-800-525-2774
Monday-Friday
7 am - 5 pm MST



Chat or E-mail
service@waterpik.com



Scan Me to
Get Started



Need more whitening tablets?

Whitening Tablet refills are available for purchase at store.waterpik.com.

Quick Start
Guide Inside!

